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INTRODUCTION

The peqGOLD MicroSpin Cycle Pure system is a convenient system for fast and reliable purification of DNA from PCR reactions within 5 minutes. The method uses PerfectBind technology to recover DNA in a range of 60 bp to 40 kb free of oligonucleotides, nucleotides polymerase and other enzymes in yields exceeding 80 %.

DNA purified with the peqGOLD MicroSpin Cycle Pure Kit is directly suitable for ligations, PCR sequencing, restriction digestion or various labelling reactions.

PRINCIPLE

Binding conditions are adjusted by addition of a specially formulated buffer, and the sample is applied to a PerfectBind DNA MicroSpin DNA column where DNA will reversibly bind to the silica matrix. Salts, free nucleotides, oligonucleotides and polymerases are removed. The DNA is eluted with relative small volume of 10 – 20 µl Elution Buffer.

Benefits

- ! Speed - DNA recovery from enzymatic reactions < 5 min.
- ! Low elution volume – Special designed column allows purification of concentrated DNA.
- ! Reliability - optimized buffers guarantee pure DNA.
- ! Safety - No organic extractions.
- ! Quality - purified DNA suitable for any application.

BINDING CAPACITY

Each MicroSpin PerfectBind DNA column can bind ~10 µg DNA.

KIT COMPONENTS

peqGOLD MicroSpin Cycle Pure Kit	5 Purifications	50 Purifications	200 Purifications
Order No.	12-6293-00	12-6293-01	12-6293-02
Components			
MicroSpin PerfectBind DNA Columns	5	50	200
2 ml Collection Tubes	5	50	200
Buffer CP	3 ml	30 ml	120 ml
DNA Elution Buffer	2 ml	15 ml	60 ml
Instruction manual	1	1	1

Buffer CP contains chaotropic salts, which are irritant. Wear gloves and other appropriate laboratory safety measures when handling.

STORAGE AND STABILITY

All peqGOLD MicroSpin Cycle-Pure Kit components are guaranteed for at least 12 months from the date of purchase when stored at 22 – 25 °C. Under cool ambient conditions crystals may form in Buffer CP. Simply warm to 37 °C to dissolve.

BEFORE STARTING

Briefly examine this booklet and become familiar with each step. Prepare all components and have the necessary materials ready before starting. peqGOLD Kits are designed to be simple, fast and reliable if all steps are followed diligently.

! All steps must be carried out at room temperature.

PEQGOLD MICROSPIN CYCLE PURE KIT PROTOCOL

Materials required, but not supplied:

! Sterile 1.5 ml centrifuge tubes

1. Agarose gel electrophoresis

Perform agarose gel/ethidium bromide gel-electrophoresis to analyze PCR product.

2. Load and Bind

Determine the volume of the PCR reaction, transfer to a clean 1.5 ml microfuge tube and add 500 μ l of Buffer CP to each 50 μ l of the PCR reaction. Vortex thoroughly.

Apply a maximum of 700 μ l of the sample to a PerfectBind DNA spin-column assembled in a clean 2 ml collection tube (provided) and centrifuge in a microcentrifuge at 10.000 x g for 2 min at room temperature. Discard the liquid. If the sample volume is more than 700 μ l, load the remaining sample to the column and centrifuge as above.

3. Elution

Place column into a clean 1.5 ml microcentrifuge tube. Add 10 – 20 μ l DNA Elution Buffer directly onto the column matrix and incubate for 1 min at room temperature. Centrifuge 1 min at 8.000 x g to elute DNA.

This represents approximately 80 – 90 % of bound DNA. An optional incubation of 5 min will increase the DNA concentration.

YIELD AND QUALITY OF DNA

Determine the absorbance of an appropriate dilution of the sample at 260 nm and at 280 nm. The DNA concentration is calculated as follows:

$$\text{DNA concentration } (\mu\text{g/ml}) = \text{Absorbance}_{260} \times 50 \times \text{Dilution Factor}$$

The ratio of $A_{260/280}$ is an indicator for nucleic acid purity. A value 1.8 – 2.0 indicates a purity of 90 % - 100 %.

Alternatively, quantity (as well as quality) can sometimes best be determined by agarose gel/ethidium bromide electrophoresis by comparison to DNA samples of known concentrations.

SHORT PROTOCOL FOR EXPERIENCED USERS

1. Determine volume of reaction. Add 500 μl volume of Buffer CP to 50 μl PCR reaction.
2. Apply solution to PerfectBind MicroSpin™ DNA column assembled in 2 ml collection tube.
3. Centrifuge at 10.000 x g for 2 min at room temperature. Discard liquid.
4. Place column into clean 1.5 ml tube and elute DNA with 20 μl DNA Elution Buffer. Centrifuge 1 min at 8.000 x g.

ORDERING INFORMATION

For extraction of DNA from agarose gels and PCR reactions:

peqGOLD Gel Extraction Kit (Classic-Line) <i>(DNA from agarose gels)</i>	12-2501-00	5 Preparations
	12-2501-01	50 Preparations
	12-2501-02	200 Preparations
<hr/>		
peqGOLD Gel Extraction Kit (Safety-Line) <i>(DNA from agarose gels)</i>	12-2500-00	5 Preparations
	12-2500-01	50 Preparations
	12-2500-02	200 Preparations
<hr/>		
peqGOLD Cycle-Pure Kit (Classic-Line) <i>(DNA from PCR reactions)</i>	12-6493-00	5 Preparations
	12-6493-01	50 Preparations
	12-6493-02	200 Preparations
<hr/>		
peqGOLD Cycle-Pure Kit (Safety-Line) <i>(DNA from PCR reactions)</i>	12-6492-00	5 Preparations
	12-6492-01	50 Preparations
	12-6492-02	200 Preparations
<hr/>		
peqGOLD MicroSpin Gel Extraction Kit (Safety-Line) <i>(DNA from agarose gels)</i>	12-6294-00	5 Preparations
	12-6294-01	50 Preparations
	12-6294-02	200 Preparations
<hr/>		
peqGOLD MicroSpin Cycle-Pure Kit (Safety-Line) <i>(DNA from PCR reactions)</i>	12-6293-00	5 Preparations
	12-6293-01	50 Preparations
	12-6293-02	200 Preparations

TROUBLESHOOTING TIPS

Problem	Likely cause	Suggestion
Low DNA yields.	Too little Buffer CP added to sample.	Add Buffer CP as indicated in the protocol.
No DNA eluted	Elution Buffer could not pass membrane	Add Elution Buffer directly to the center of the column.
	Too little Buffer CP added to sample.	Add Buffer CP as indicated in the protocol.
Optical densities do not agree with DNA yield on agarose gel.	Trace contaminants eluted from column increase A_{260}	Add more Buffer CP to PCR reaction mix.

D PEQLAB Biotechnologie GmbH, 91052 Erlangen, Freecall (D): 0800 100 20 16, info@peqlab.de, www.peqlab.de
AT PEQLAB Biotechnologie GmbH, 6404 Polling, Tel: +43 (0) 5238 84 169, info@peqlab.at, www.peqlab.at
UK PEQLAB Ltd., Fareham PO15 5TT, Freecall (UK): 0808 202 1302, info@peqlab.co.uk, www.peqlab.co.uk
USA PEQLAB LLC, Wilmington, DE 19810, Toll-Free (US): 877 737 5220, info@peqlab.us, www.peqlab.us

Creating the future together.