



Certified Water



Water is Life.....

MO BIO is committed to helping people better their lives by contributing resources to provide clean, accessible water. MO BIO donates a percentage of the sales of all our water products to Water Partners International. For more info go to www.water.org



This symbol indicates that all products on this page are part of our Water is Life Program. (See page 78)

PCR Water	
Catalog No.	Quantity
17000-1	1 ml
17000-5	5 x 1 ml
17000-10	10 x 1 ml
17000-11	10 ml bottle
Molecular Biology Grade Water	
Catalog No.	Quantity
17012-200	200 ml
17012-5200	5 x 200 ml
DEPC Treated Water	
Catalog No.	Quantity
17011-200	200 ml
17011-5200	5 x 200 ml
Endotoxin-Free Water	
Catalog No.	Quantity
17013-10	10 ml
17013-50	50 ml
17013-100	100 ml
17013-500	500 ml

PCR Water

Water certified by PCR to be free of all DNA

DESCRIPTION PCR Water serves as a reliable contaminant-free source for PCR experiments. Manufactured using a proprietary procedure that removes all detectable DNA contamination, PCR Water is certified to be free of all DNA, RNase, and DNase contamination. In addition, it is validated for use with highly sensitive PCR reactions including those with eubacterial universal primers.

Features and Benefits

- Made specifically for use as a component for PCR reactions
- Double pouched for additional contaminant protection
- Sterile and UV irradiated
- Validated for use with highly sensitive PCR reactions

Molecular Biology Grade Water

High quality water that is certified RNase and DNase free

DESCRIPTION Certified to be free of all RNase and DNase contamination. This water serves as a reliable contaminant-free source for DNA and RNA experiments.

DEPC Treated Water

DEPC treated water that is certified RNase and DNase free

DESCRIPTION DEPC Treated Water is produced specifically for sensitive laboratory applications, is sterile and certified to be free of all RNase and DNase contamination. Treated with DEPC, this water serves as a reliable contaminant-free source for experiments.

Endotoxin-Free Water

Water that is certified free of endotoxins

DESCRIPTION Endotoxin-Free Water is produced specifically for sensitive applications where low endotoxin levels are critical. This water source is ideal for preparing solutions and reagents for transfection and injection procedures and is certified to contain <0.001 EU/ml of water. Endotoxin-Free Water is also useful for elution of endotoxin-free plasmid DNA.