

SHAKER AND STIRRER



V1 PLUS

BOECO BOECO VORTEX MIXER V 1 PLUS

Vortex mixer V1 plus with variable speed control. Continuous operation as well as with automatic starter, triggered by pressing the attachment. Incl. standard attachment for test tubes up to 20 mm Ø

Specification

| | |
|-----------------------|---------------------------------------|
| Kind of motion: | orbital, 4 mm |
| Speed range (1/min): | 300 - 3000 1/min |
| Time of acceleration: | 3 s |
| Tube's capacity: | 1,5-50ml |
| Speed display: | scale (min/max) |
| Switch: | touch / continuous mode |
| Overall size: | 80 x 140 x 70 mm |
| Weight: | 0,9 kg |
| Power: | 12 V. 500 mA (External power adapter) |

| Code | Description |
|-------------|---------------------------------------|
| BOE 8055000 | Vortex mixer V1 plus, AC Adapter 220V |
| BOE 8055100 | Vortex mixer V1 plus, AC Adapter 110V |



MMS 3000



BOECO MAGNETIC STIRRER MMS 3000

Motor drive magnetic stirrer MMS-3000, with stainless steel plate and speed regulation up to 3000 rpm. A optional to be ordered detachable support rod allows to fix different sensors (temperature, pH etc. inside the Liquid).

Specification

| | |
|-----------------------|---------------------------------------|
| Speed range: | 0 - 3000 1/min |
| Max. stirring volume: | up to 2,0 L |
| Platform size: | 150 mm diam. |
| Platform material: | stainless steel |
| Support rod height: | 320 mm |
| Overall size: | 165 x 200 x 65 mm |
| Weight: | 1,2 kg |
| Power: | 12 V. 300 mA (External power adapter) |

| Code | Description |
|-------------|--|
| BOE 8056500 | Magnetic stirrer MMS-3000, AC Adapter 220V |
| BOE 8056510 | Magnetic stirrer MMS-3000, AC Adapter 110V |
| BOE 8056501 | Detachable support rod |



MSH 300

BOECO MAGNETIC STIRRER WITH HOTPLATE MSH 300

Motor drive Magnetic stirrer with aluminum hot plate for simultaneously stirring and heating of chemical reagents. Equipped with feed-back thermoregulation.

Specification

| | |
|-------------------------|-------------------|
| Speed range: | 250 - 1250 1/min |
| Temperature regulation: | 30 - 330 C°. |
| Max. stirring volume: | Up to 2,0 L |
| Hot Platform size: | 150 x 150 mm |
| Hot Platform material: | aluminium |
| Overall size: | 170 x 210 x 95 mm |
| Heating power: | 600 Watts |
| Weight: | 2,5 kg |

| Code | Description |
|-------------|--|
| BOE 8057000 | Magnetic stirrer with hotplate MSH 300, 220V |
| BOE 8057100 | Magnetic stirrer with hotplate MSH 300, 110V |