



MR-1

## BOECO MINI-ROCKER SHAKER MR-1

provides gentle rocking motion. It is perfect for soft mixing of biological components. Incl. black non-slip rubber mat

### Specification

Kind of motion: rocking with fixed 7° angle  
 Speed range: 7,5 - 30 rpm  
 Speed display: scale (min/max)  
 Timer: digital, 1 min - 23h. 59 min.  
 Platform size: 200 x 200 mm  
 Overall size: 220 x 205 x 120 mm  
 Max Loading: 500 g  
 Weight: 1,6 kg  
 Power: 12 V. 500 mA (External power adapter)

Code	Description
------	-------------

BOE 8058500	MR-1, Mini-rocker Shaker, AC Adapter 220V
BOE 8058510	MR-1, Mini-rocker Shaker, AC Adapter 110V

PSU-2T  
PLUS

## BOECO MICROTITER PLATE SHAKER PSU-2T PLUS

specially designed for medical diagnostic. Orbital shaking platform IPP-2 with a holding mechanism for 2 standard microtiter plates. The variable speed control provides a optimally adapted speed range from 150 to 1000 rpm with a orbit of 2 mm. .

### Specification

Orbit size: 2 mm  
 Speed range: 150 - 1000 1/min  
 Speed display: scale (min/max)  
 Timer: digital, 1 min - 23h. 59 min.  
 Overall size: 220 x 205 x 90 mm  
 Weight: 2,5 kg  
 Power: 12 V. 500 mA (External power adapter)

Code	Description
------	-------------

BOE 8052000	Microtiter Plate Shaker PSU 2T plus, AC adapter 220V
BOE 8052100	Microtiter Plate Shaker PSU 2T plus, AC adapter 110V
BOE 8052002	Platform IPP-4 for 4 Microtiter Plates

PST-60 HL  
PLUS

## BOECO MICROPLATE INCUBATOR SHAKER PST-60 HL PLUS

The PST-60 HL plus provides orbital shaking of the samples and has a heating platform for 2 x 96 wells microplates and an upper heating lid. The display indicates simultaneous set and actual time, temperature and shaking speed. The PST-60 HL plus is ideal for a variety of cytochemistry, immunochemistry, biochemistry and research applications.

### Specification

Temperature regulation range: 25°C ..... 60 °C  
 Accuracy: ± 0,1 °C  
 Temp. Uniformity over the platform: ± 0,2 °C  
 Time of thermoblock heating from RT till 37°C: 15 - 20 min.  
 Orbit : 2 mm  
 Speed regulation: 250 -1200 rpm / increment 10 rpm  
 Timer with sound signal: 0 min - 96 hrs, increment 1 min  
 Platform dimensions: 250 x 150 mm  
 Platform capacity: 2 x 96 wells microplates  
 Power (external) : DC 12 V. 4,16 A  
 Size : 270 x 260 x 125 mm  
 Weight : 6,0 kg

Code	Description
------	-------------

BOE 8072000	Microplate Incubator Shaker PST-60 HL plus, AC 100-240V 50/60 Hz
-------------	--