



OS-20 WITH PP-4 PLATFORM

## BOECO UNIVERSAL ORBITAL SHAKER OS-20

Compact, low-profile, orbital laboratory shaker for versatile applications in the fields of biotechnology, microbiology and medical diagnostics. The OS-20 provides gentle rotational mixing and precise circulation at each point of the platform. Variable digitally displayed speed control from 50-250 rpm. Digital timer for setting a time interval and continuous operation.

### Specification

Kind of motion: orbital  
 Orbit size: 20 mm  
 Speed regulation: 50 - 250 rpm  
 speed display: digital (LED)  
 Timer: from 1 min to 999 min / continuous  
 Time display: digital (LED)  
 Max load: 2,5 kg excl. attachment  
 Power: 12V. 500 mA (External power adapter)  
 Overall size  
 (w/o attachment): 265 x 270 x 70 mm  
 weight: 3,9 kg.

Code	Description
BOE 8059000	Universal Orbital Shaker OS-20, AC Adapter 220V
BOE 8059100	Universal Orbital Shaker OS-20, AC Adapter 110V



OS-10 WITH UP-12 PLATFORM

## BOECO UNIVERSAL ORBITAL SHAKER OS-10

compact, powerful and quite microprocessor controlled orbital laboratory shaker for many general purpose shaking application.


Speed and time setting are easy and during operation the LCD status display clearly indicates both set and actual speed, plus set and elapsed time.

### Specification

Kind of motion: orbital  
 Orbit size: 10 mm  
 Speed regulation: 50 - 350 rpm  
 speed display: digital (LCD)  
 Timer: from 1 min to 96 hours / continuous  
 Time display: digital (LCD)  
 Max load: 3 kg excl. attachment  
 Power: 12V. 500 mA (External power adapter)  
 Overall size  
 (w/o attachment): 270 x 260 x 85 mm  
 weight: 3 kg

Code	Description
BOE 8069000	Universal Orbital Shaker OS-10, AC Adapter 220V
BOE 8069100	Universal Orbital Shaker OS-10, AC Adapter 110V

## ACCESSORIES FOR OS-20 AND OS-10



PLATFORM P-6/250

Code	Description
BOE 8059001	Dish platform PP-4 with rubber mat for petri dishes, culture flasks etc.
BOE 8059002	Universal platform UP-12 for all kind of flasks 260 x 180 mm, stainless steel with variable clamping rolls and rubber mat. Numbers of flasks at the platform 100- 150 ml flasks max 12-9 pieces 200- 250 ml flasks max 6-4 pieces 500-1000 ml flasks max 4-2 pieces
BOE 8059003	Platform P-6/250 with fixed clamps for 6 flasks 250/300 ml.