


Cryo Products Cryo Tubes



Cryo Tubes

Polypropylene with highest translucence.
For the storing and conservation of tissue cultures and biological samples e.g. sera, blood and spermata.
For applications to -190 °C in the gaseous phase of liquid nitrogen or in freezers.

- pyrogene-free
- DNase-, RNase- and DNA-free
- conical bottom inside for optimum sample recovery
- white label and white graduation
- Cryo Tubes with male thread
- positive leakproof screw cap available in two version:
 - with integrated lips without additional sealing ring
 - with silicone-o-ring sealing
- Cryo Tubes with self standing base provide easy one-hand operation in the Work Stations
- 2 ml Cryo Tubes also available with round base

Packed in cardboard boxes with 10 PE-bags x 100 pcs.

Volume	Description	Size Dia. x Hight	Quantity per Pack	Order Number
Cryo Tubes with integrated Lip Sealing				
1.2 ml	self standing	12.5 x 43 mm	1000	51 10 012
2.0 ml	round base	12.5 x 48 mm	1000	51 10 020
2.0 ml	self standing	12.5 x 49 mm	1000	51 10 022
3.0 ml	self standing	12.5 x 72 mm	1000	51 10 030
4.0 ml	self standing	12.5 x 78 mm	1000	51 10 040
5.0 ml	self standing	12.5 x 92 mm	1000	51 10 050
Cryo Tubes with Silicone O-Ring Sealing				
1.2 ml	self standing	12.5 x 43 mm	1000	51 10 112
2.0 ml	round base	12.5 x 48 mm	1000	51 10 120
2.0 ml	self standing	12.5 x 49 mm	1000	51 10 122
3.0 ml	self standing	12.5 x 72 mm	1000	51 10 130
4.0 ml	self standing	12.5 x 78 mm	1000	51 10 140
5.0 ml	self standing	12.5 x 92 mm	1000	51 10 150