

2.6 Gelloader Tips



1 - 200 ul Volume Range

Free of human DNA, certified free of DNase, RNase and PCR inhibitors



- The new and unique AHN capillary tips
- Fine point gel tips for extreme precision, Gel 0.7 mm
- These round tips are very useful for SDS-PAGE gels
- Universal fitting for most pipette brands
- Injection moulded from virgin polypropylene in low retention quality
- The new and unique AHN capillary filter tips
- Useful for when one set of pipettes is used for both experimental set-up and results
- Universal fitting for most pipette brands
- The most effective filter for tips. Pore size effectivity 4-6 microns.

Specification	According to
General requirements	ISO Standard 8655
CE-marking certificate	EU-Directive 98/79/EG, competent authority code DE/CA90
Sterility certificate	EN 552/ISO 11137 and DIN EN 556
RoHS Compliance certificate	EU Directive 2002/95/EC
Environmental Management certificate	ISO Standard 14000 and ISO Standard 14001
USP-Plastic Test, Class VI, certificate	ISO Standard 10993
Conformity testing certificate	DIN Standard 12600
GMP/GLP certificate	WHO Technical Report Series 908, 37. Report, Geneva 2003

Cat.no.	Description	Sterile	Colour	Pack type	Pack size	Sales unit
1-130-50-0	AHN Gelloader Tips 1-200 ul	non	natural	bag	1000	1000
1-131-50-0	AHN Gelloader Tips 1-200 ul	yes	natural	bag	1000	1000
1-130-96-0	AHN Gelloader Tips 1-200 ul	non	natural	rackbox	96	10 x 96
1-131-96-0	AHN Gelloader Tips 1-200 ul	yes	natural	rackbox	96	10 x 96
2-130-96-0	AHN Gelloader Filtertip 30 ul	non	natural	rackbox	96	10 x 96
2-131-96-0	AHN Gelloader Filtertip 30 ul	yes	natural	rackbox	96	10 x 96
2-132-96-0	AHN Gelloader Filtertip 100 ul	non	natural	rackbox	96	10 x 96
2-133-96-0	AHN Gelloader Filtertip 100 ul	yes	natural	rackbox	96	10 x 96